

Fluorescence immunostaining of vibratome sections

Original protocol from Huttner Lab, modified by Zerial lab

- 1) Dissect the liver and cut into pieces (~5x5x5 mm)
- 2) Fix the liver pieces with 4% PFA in PBS for 2 hours at room temperature then overnight at 4°C
- 3) Wash liver pieces in PBS
- 4) Embed the liver lobes into 4% agarose in PBS
- 5) Make vibratome sections (50 ~ 200 µm thickness)
- 6) Place selected slices in 48-welled plate (1 slice per well).
- 7) Permeabilize with 0.5% Triton X-100 in PBS for 60 min at room temperature.
- 8) Quench sections in 10mM NH₄Cl in PBS for 30 min at room temperature.
- 9) Wash with TxBuffer 3 times for 3 min each.
- 10) Add primary antibody in TxBuffer (2 overnights at 4°C).
- 11) Wash with 0.3% Triton/PBS 5x 5 min.
- 12) Add secondary antibody in TxBuffer (2 overnights at 4°C).
- 13) Wash with 0.3% Triton/PBS 5x 5 min.
- 14) Wash in PBS 3x 1 min.
- 15) Clear the samples with glycerol (to reduce background):
 - Incubate in 25% glycerol/PBS until samples settling on the bottom of dish
 - Increase the glycerol concentration to 90% in a step-wise manner (50% → 75% → 90%):
 - i. For overnight change to 50% glycerol in PBS
 - ii. Next morning change to 75% glycerol in PBS
 - iii. For overnight change to 90% glycerol in PBS
- 16) Mount on glass slide with 90% glycerol/0.2% n-propyl-gallate covered by #1.5 coverslips (thickness 0.17 ± 0.005 mm)
- 17) Image

TxBuffer (1L):

0.2% gelatin (100mL 2% gelatin in PBS filtered, stored at -20°C)

300mM NaCl (17.53g)

0.3% Triton X-100 (3mL)

Add PBS to make 1L, aliquot to 50mL falcon tubes, store at -20°C