

### Isolation of primary hepatoblasts

Protocol adapted from Tanimizu N et al., *J Cell Sci.* 2003 and Kamiya A et al, *EMBO J* 1999

- 1) Remove the livers from E14.5 embryos
- **Digestion of the embryonic livers:**
- 2) Place embryonic livers into prewarmed Liver Perfusion Medium (LPM) (Invitrogen, No.17701) (1.5 mL medium per liver).
- 3) Break livers mechanically by pipeting the tissue in LPM 3-5 times through a glass pipette.
- 4) Incubate livers in LPM for 20 min at 37 °C waterbath, mix with turning after every 5 min.
- 5) Centrifuge at 212×g, 3 min at room temperature.
- 6) Remove supernatant and loosen the pellet (necessary for avoiding aggregation).
- 7) Suspended liver pieces in 12ml Liver Digest Medium (LDM) (Invitrogen, No.17703-034) supplemented with DNase I (1:1000).
- 8) Incubate liver pieces in LDM for 20 minutes at 37 °C, pipetting every 5 minutes.
- 9) Drain cell suspension through a cell strainer (70 μm) (BD Falcon, No.352350). Add PBS until 40 mL of cell suspension.
- 10) Centrifuge at 260×g, 3 min at room temperature.
- 11) Remove the supernatant and loosen pellet to avoid aggregation.
- 12) Resuspend cells with 10 mL of PBS.
- 13) Drain the suspension through cell strainer (70 μm) again and wash the tube and the cell strainer with additional 30 mL of PBS.
- **Hemolysis**
- 14) Centrifuge at 260×g for 3 min at room temperature.
- 15) Remove supernatant and loosen pellet.
- 16) Add 12 mL of ice-cold Hemolysis buffer (12 mL per 8 livers) and resuspend the pellet.
- 17) Incubate on ice for 5 min.
- 18) Gently add 23 mL ice-cold Hepatoblasts Differentiation Medium (HDM) containing 5% FBS. Remove debris by passing suspension through 70 μm cell strainer.
- 19) Wash the walls of the tube with additional 5 mL of ice-cold HDM and pass the solution through the cell strainer.
- 20) Centrifuge at 260×g for 3 min at room temperature.
- 21) Remove supernatant and loosen pellet.
- 22) Resuspend cells with 10 mL ice-cold HDM and pass through the cell strainer.
- 23) Wash the tube again with 30 mL PBS and drain through cell strainer.
- 24) Count the number of cells.
- 25) Centrifuge at 260×g for 3 min at room temperature.
- 26) Remove supernatant and loosen pellet.
- 27) Resuspend cells (100,000,000 cells / mL) into MACS buffer and split the cell suspension between eppendorfs (1.5 mL), 300 μL of cells in MACS buffer into one eppendorf.
- **Magnetic cell sorting**
- 28) Incubate cells with anti-FcγR antibody (1:100 dilution) on ice for 10 min.
- 29) Add FITC anti-Dlk1 antibody (1:40 dilution) and incubate on ice for 15 min.
- 30) Add 1 mL MACS buffer and centrifuge at 660×g (tabletop centrifuge) for 3 min at 4°C.
- 31) Remove supernatant and loosen pellet.
- 32) Wash again with 1 mL MACS buffer and centrifuge at 660×g (tabletop centrifuge) for 3 min at 4°C .

- 33) Remove supernatant and loosen pellet.
- 34) Resuspend cells with 150  $\mu$ l MACS buffer / 30,000,000 cells and add anti-FITC-microbeads (1:10 dilution). Incubate on ice for 15 minutes.
- 35) Add 1 mL MACS buffer and centrifuge at 660 $\times$ g (tabletop centrifuge) for 3 min at 4°C.
- 36) Remove supernatant and loosen pellet.
- 37) Resuspend cells in 3 mL MACS buffer. Remove debris by passing through a cell strainer (70  $\mu$ m) into a 15 mL tube.
- 38) Centrifuge at 260 $\times$ g for 3 min at room temperature.
- 39) Remove supernatant and loosen pellet.
- 40) Resuspend in 500  $\mu$ l of MACS buffer (up to 100,000,000 cells can be passed through the MACS column).
- 41) Install a MACS column (MS column, Miltenyi Biotech, No.130-042-201) onto a magnetic holder and equilibrate it by washing through 500  $\mu$ l MACS buffer.
- 42) Load the 500  $\mu$ l cell suspension onto the MACS column and collect flow.
- 43) Load again.
- 44) Wash column 3 times with 500  $\mu$ l MACS buffer.
- 45) Take MACS column off magnetic holder.
- 46) Put MACS column onto a 15 mL tube and elute Dlk1-positive cells with 3x 1 mL MACS buffer.
- 47) Count the number of cells.
- 48) Centrifuge at 260 $\times$ g for 3 min at room temperature.
- 49) Remove supernatant and loosen pellet.
- 50) Resuspend cells in icecold HDM (1,000,000 cells / mL)

Hepatocyte differentiation medium:

- DMEM (High glucose 4.5 g/mL)
- 1x Pen/Str
- 1x L-glutamine (Invitrogen, No. 25030-024)
- 1x non-essential amino acids (Invitrogen, No. 11140-035)
- Dexamethasone 31.25  $\mu$ L per 500 mL of medium (Sigma-Aldrich, D1756-25MG)
- Remove 25 mL of sf medium
- add 25 mL of filtered serum (this serum is not heat-inactivated!) (filter and then immediately wash the filter with some medium to get all serum proteins)

MACS buffer:

- 0.5% BSA
- 2 mM EDTA
- PBS
- filter the buffer before use!

